



Glass Oven G-300

The Multipurpose Solution For Small Samples

Only \$11,250

Usually \$12,500

Code: 11G3001

Angle adjustment button


To adjust the angle of the heating tube.

SET button

To set heating temperature.

Heating tube

With semiconductor coating for gentle heating and fast heating and cooling times.



STOP button

To stop the process and the connected vacuum pump.

Interface

Showing heating temperature and rotation speed.

Navigation control

To adjust heating temperature and rotation speed.

Heating temperature:

up to 300 °C

Volume:

Drying: up to 250 mL (sample vol.),
Distillation: 5–40 mL (flask size)

Introducing the Glass Oven G-300, BUCHI's all-in-one instrument for small samples.

With interchangeable accessories, the Glass Oven G-300 can perform distillation, drying, sublimation, and freeze-drying, eliminating the need for multiple large instruments.

Its compact size even allows placement in a glovebox, ensuring the safe handling of air-sensitive samples.



Ref. Code: BUCHI-G-300

Valid Until: 31 Aug 2025



1800 010 114
customerservice.ins@dksh.com
labshop.dksh.com.au/shop/



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.