



Printed Culture Tube

Only \$40

Usually \$113

KIMAX® culture tube has a fire-polished top rim and uniform wall thickness for maximum heat transfer and chemical resistance.

Code:
KIML45048-16150_U
72 items

- Complies with ASTM Specification E982, Type IV requirements
- Manufactured from borosilicate glass 3.3
- Durable white ceramic enamel marking spot

Plain Borosilicate Glass Tubes



Code	Size		Usually	Only
KIML73500-16150	23ml (16mm x 150mm)	1000 items	\$316	\$80
KIML73500-650	1ml (6mm x 50mm)	1000 items	\$178	\$30
KIML73500-16125	19ml (16mm x 125mm)	1000 items	\$293	\$161
KIML73500-18150	28.5ml (18mm x 150mm)	500 items	\$236	\$180

These premium quality tubes feature sturdy, uniform bottoms and consistent lengths.

- Complies with ASTM Specification E890
- 51 expansion borosilicate glass offers excellent chemical resistance against sodium leaching, a factor common in soda lime glass products
- Tubes are packed in convenient trays for ease of use on bench tops and in drawers
- Trays are shrink-wrapped in plastic film, with tubes remaining in constant alignment for protection against in-transit breakage and contamination
- No marking spot

Special pricing valid for instock items only

Ref. Code: DWK-Tubes

Valid Until: 31 Jul 2025