



79 Litre Top-loading Autoclave

Designed To Fit Your Space And Needs

Only \$10,975

Usually \$17,342

Code: TOMYSX-700E

Including two baskets




TOMY SX-Series Autoclaves are floor-standing lab autoclave sterilizers, a design that saves you valuable bench-top space.

The slim profile creates a small footprint that maximizes your floor space, while the large capacity fits all types of lab media, glassware, liquids, agar, etc.

This compact, high-capacity lab autoclave helps you make more room for research.

- Simple lid opening using either hand or foot
- Safety devices include: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection & safety valve
- Warming cycle enabling media to be kept in liquid state for dispensing
- Maximum sterilising temperature to 135 °C
- Removable condensate bottle, negating the need for drains
- 2-year warranty



- Chamber dimensions:
370mm Ø x 774mm H
- Max Bottle Capacity:
24 x 1L, 33 x 500 mL, 51 x 250 mL
- Max Flask Capacity:
12 x 1L, 21 x 500 mL, 33 x 250 mL

Ref. Code: TOMY-SX700E

Valid Until: 31 Jul 2025



1800 010 114
customerservice.ins@dksh.com
labshop.dksh.com.au/shop/



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.