



designed for scientists

C-Mag HP 4 Hotplate

In Stock

Only \$207

Usually \$1,125

Code: IKAA3581605

**Only 6 units
available at this
super low price!**



Hotplate made of ceramic glass offering excellent chemical resistance.

- Fixed safety circuit of 550 °C
- Hot Top indicator >> hot surface warning to prevent burns!
- Exact temperature setting via digital display (LED)
- Digital error code display
- Elevated control panel for protection against leaking liquids

Ref. Code: IKAA-CMAG-HP4

Valid Until: 30 Sep 2025



1800 00 84 53
sales.au@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.com.au/shop/



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.