



Miele

PG 8583 ADP Washer/Disinfector

Only \$12,126

Code: MIEL62858306

Included **FREE**

1 x A102 Upper Basket with spray arm
1 x A150 Lower Basket for optional modules


Valued at \$1,517



A102 Upper Basket



A150 Lower Basket



PLUS 15% discount on optional accessories with free Detergent & Neutraliser – please call us to discuss your complete washing solution

- Powerful 500L/min pump with integrated heating elements
- In-built steam condenser
- Spray arm & spray pressure monitoring ensuring a clean wash
- Capable of washing up to 128 flasks or 98 pipettes with optional accessories
- 14 programs and 2 free program positions
- Multi-stage filtration system
- Stainless steel interior/exterior for durability & cleanliness
- Passive EcoDry drying function – at 70°C door opens slightly enabling the load to be dried
- ADP model including pump for AD water
- Two wash levels
- Dispensing pumps for liquid detergent and liquid neutraliser ensure accurate dosing
- 3 phase power, convertible to single phase (charges apply)
- 2-year warranty

Ref. Code: MIEL-PG8583

Valid Until: 31 Aug 2025

 **bio-strategy**
Part of DKSH Group



1800 00 84 53
sales.au@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.com.au/shop/

 **DKSH**

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.