

No expensive drainage required!

Integrated condensate bottle negates the need for connection to drains.




**Only \$13,995**

**Usually \$20,965**

Code: TOMYFLS-1000

Includes:

2 stainless steel mesh baskets  
& 2-year warranty



Easy opening foot  
pedal operated lid

- Full safety features including water level detection and lid interlock
- Large 100L vertical chamber capacity
  - » Max Bottle Capacity: 34 (17 x 500mL per basket)
  - » Max Flask Capacity: 4 x 5,000mL
- Chamber dimensions: 450 Ø x 678 mm H
- External dimensions (mm): 680 (D) x 680 (W) x 1,022 (H) (including the pedal and control panel)
- Maximum sterilising temperature to 123 °C
- Weight: 108kg

Good range of accessories available including:  
punched metal stainless steel baskets, solid stainless steel baskets,  
printer, external sensor (wandering probe) and preventative maintenance contracts.

Ref. Code: TOMY-FLS

Valid Until: 31 Aug 2025