



Melting Point M-565

Melting Point At A Touch Of A Button

Only \$28,000

Usually \$31,000

Code:11058004M-565



The video camera equipped Melting Point M-565 for the automatic determination of melting and boiling points is designed for true, professional quality control.

When used in combination with the Sample Loader M-569, determinations are achieved with unmatched reproducibility.

Features:

- Fully automatic melting point determination
- Real time video displayed in 6 x magnification
- Video replay function with variable speed
- Temperature range between room temperature and 400 °C
- Large lens with 2.5 magnification for precise determination of the samples
- Heating block protective coating allows for easy cleaning and extended lifetime
- Compliant with Pharmacopeia methods including PH. EUR. 6.1.2.2.60, USP XXI 741 and JP
- GLP/GMP compliant protocols either with printer

Ref. Code: BUCHI-M565

Valid Until:31 Oct 2025



1800 010 114
customerservice.ins@dksh.com
labshop.dksh.com.au/shop/



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.