




# Superfrost Plus White Slides

POWERED BY SHANDON MICROM MENZEL-GLÄSER Richard-Allan Scientific

Superfrost PLUS slides designed to secure tissue sections without the need for laboratory-applied coatings

**Only \$32**

Code: EPBRSF41296SP



Electrostatically charged adhesion slides with excellent adhesive properties for high reliability in instrument-based applications.

- Uniform coating
- Minimal sample loss
- For tissue slices 2-5µm thick
- Eliminates background staining in standard H&E stains

**Compatible with:** Epredia Microtomes and Cryostats



Contact us for pricing

Ref. Code: EPBR-SFSPWH

Valid Until: 31 Aug 2025



1800 00 84 53  
[sales.au@bio-strategy.dksh.com](mailto:sales.au@bio-strategy.dksh.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
[labshop.dksh.com.au/shop/](http://labshop.dksh.com.au/shop/)



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.