

**Limited  
Stock  
Available**

ROCKER 2D digital

**Only \$1,472**

Usually \$2,330

Code: IKAA4003000



Anti-slip rubber mat included

Digital rocking shaker with adjustable speed for smooth mixing tasks.

Used for homogeneous mixing in flasks, culture flasks, Petri dishes and tubes.

Even with load and voltage fluctuations, shaking speed remains constant.

An adjustable tilt angle ensures optimum adaption to the particular experimental conditions.

- Timer: Countdown from 1s to 99h 59min
- Counter: Display of shaking time
- Easy operation with touch keypad
- Shaking diameter: 0 - 15° (angle)
- Speed range: 5 - 80 rpm
- Error Code Display
- Usable at 4 to 50 °C in incubators
- Suitable for continuous operation
- Different attachments enable use for a variety of applications

KS 501 Digital  
Orbital Shaker

**Calls us for a  
competitive price**

Code: IKAA4003000



Low profile laboratory orbital shaker with large mounting surface and load capacity of up to 15 kg.

Ideal for biological and microbiological growth tests and for production of stock solutions.

Suitable for vessels with a volume of more than 250 ml, e.g. round flasks, Erlenmeyer flasks, culture flasks and bottles.

- Infinitely variable speed control from 0-300 rpm
- Digital display
- Orbital diameter 30 mm
- Continuous operation (endless) even under extreme loads
- Timer function
- Integrated USB and RS 232 interface
- Platform size : 450 x 450 mm
- Attachments are not included, please order separately

Ref. Code: IKA-Digital

Valid Until: 30 Sep 2025