

# BASE

## Front Loading Autoclave



**Only \$43,000**

**Usually \$45,934**

Code: PRIOEH150-BASE

With cylindrical 150L chamber & 2 non-tilting removable shelves. Chamber carries a 10-year warranty

**Export Crate**

Code PRIOPC/AEC/004 **Only \$1,910** | Usually \$2,355

TACTROL 3 microprocessor control system with simple and fully variable setting of process time, temperature and with graphic indication of cycle status.

Many optional accessories available including loading trolleys, discard containers, printer, wandering probe, water softener, etc – please contact Bio-Strategy Part of DKSH to assist in configuring your autoclave solution.

- Quickseal single-action door fitted with thermal and pressure locks to prevent opening at load temperature above 80 °C and pressures above 0.2 bar
- Forced air cooling as standard helping to reduce cycle times
- Automatic timed & pulsed free steaming as standard for improved air removal and sterilising performance
- Steam generated within the chamber
- Delayed start and media warming function as standard for ready to pour media at the start of the day
- Low water level sensor
- Epoxy coated panels and frame members are treated with an anti-bacterial agent which is effective against all bacteria and fungi including MRSA
- Pressure vessel carries a 10-year warranty and is insurance approved

Ref. Code: PRIO-EH150

Valid Until: 31 Oc 2025