

Versatile Climate Chamber



Only \$14,715

Usually \$20,150

Code: SANYMLR-352H-PE

294L Capacity

Includes adjustable shelves:

4 wire & 1 solid

The MLR Climate Chamber, with LED lighting is suitable for a wide range of applications, including plant growth and insect studies. The wide variety of temperatures, humidity and lighting patterns that are essential in research, environmental studies, and testing can be accurately reproduced and controlled.

Designed with a slim cabinet profile makes for efficient use of available floor space. Robust cabinet construction with corrosion-resistant stainless steel surfaces simplifies cleaning and eliminates rust due to high humidity.

Lamps

High-efficiency LED lamps are arranged in programmable banks of five lamps in the door, left and right interior sides, and are programmable in up to 6 combinations from zero to 15 lamps on.

Programmable Setting

The temperature and light settings are programmable to provide the perfect cycles for your applications. The unit allows the programming of starting day and time of operation and 10 programmes of up to 12 steps each to be memorized. Switch between the selectable clock mode or timer mode and use of the 'Join Mode' for allowing multiple programs to be linked.



Wide Range of Applications:

The PHCbi climatic test chambers support a wide range of applications. Creating the optimum environments for research and testing protocols.

Ref. Code: SANY-MLR352

Valid Until: 31 Dec 2025