



Miele

PG8583 Washer/Disinfector 145L Capacity Under bench or Stand-alone

Clearance Price
Only \$9,700

Usually \$11,844
Includes free Detergent & Neutraliser
Code: MIEL62858304



With passive EcoDry drying function – activates as the temperature descends to 70°C, opening the door slightly, enabling the load to dry.

- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- In-built steam condenser
- Spray arm and spray pressure monitoring, ensuring a clean wash
- Capable of washing up to 128 flasks or 98 pipettes with optional accessories
- Multi-stage filtration system
- Stainless steel interior/exterior for durability and cleanliness
- Compact: 600mm x 600mm x 850mm [W x D x H] (820mm with lid removed)
- 3 phase power, convertible to single phase (charges apply)
- 2-year Warranty

PLUS 15% Off
optional accessories

**Call us to discuss your complete
washing solution**



Valid Until: 28 Feb 2026

Ref. Code: MIEL-PG8583AD

1800 008 453 | sales.au@bio-strategy.dksh.com | www.bio-strategy.com | labshop.dksh.com.au

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Delivering Growth

 **bio-strategy**
Part of DKSH Group

 **DKSH**