

# Miele

## PG8583 ADP Washer/Disinfector

**Stock Clearance**

**Further Price  
Reduction**



**Only \$9,850**

Usually \$12,126

Code: MIEL62858306



**FREE**

Detergent & Neutraliser



**Plus 15%**

discount on optional accessories

**Please call us to discuss your complete  
washing solution**

Capable of washing up to 128 flasks or  
98 pipettes with optional accessories

- Underbench or stand-alone model, 145 litre capacity
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- In-built steam condenser
- Spray arm and spray pressure monitoring ensuring a clean wash
- 14 programs and 2 free program positions
- Multi-stage filtration system
- AO value control
- Stainless steel interior/exterior for durability and cleanliness

- Compact: 600mm x 600mm x 850mm W x D x H
- Passive EcoDry drying function – activates as temperature decends to 70°C opening the door slightly enabling the load to be dried.
- ADP model including pump for AD water
- Two wash levels
- Many options and accessories to select from
- Dispensing pumps for liquid detergent and liquid neutraliser ensure accurate dosing
- 3 phase power, convertible to single phase (charges apply)
- 2-year warranty

Valid Until: 28 Feb 2026

Ref. Code: MIEL-PG8583ADP

1800 008 453 | [sales.au@bio-strategy.dksh.com](mailto:sales.au@bio-strategy.dksh.com) | [www.bio-strategy.com](http://www.bio-strategy.com) | [labshop.dksh.com.au](http://labshop.dksh.com.au)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Delivering Growth

 **bio-strategy**  
Part of DKSH Group

 **DKSH**