



# Top-loading Autoclaves 79L & 100L

**In Stock - At 2025 pricing!**

Take advantage before new 2026 stock arrives

**Only \$9,975**

Usually \$15,840

Includes 2  
punched stainless steel baskets

Code: TOMYSX-700E



Has a removable condensate bottle negating the need for drains and is powered from a 240V 50Hz 15amp supply .

- Sterilization, warming, and heating cycles
- Warming cycle enabling media to be kept in liquid state for dispensing
- Maximum sterilising temperature to 135 °C
- Dual-fan cooling system for faster cycle times
- Easy opening foot operated lid
- 2-year warranty

**Only \$13,500**

Usually \$24,000

Includes 2  
punched stainless steel baskets

Code: TOMYFLS-1000



With an integrated condensate bottle negates the need for connection to your drainage network!

- 3-level cooling fan for faster run times
- Excellent basket loading height: 785mm
- Large chamber for bulk or large volume items: 450mm Ø x 678mm H
- Sterilising temperature range from 106°C to 123 °C
- Easy-opening foot pedal operated lid
- 2-year Warranty

Good range of accessories available, solid stainless steel baskets, printer, external sensor (wandering probe) and preventative maintenance contracts.

Valid Until: 28 Feb 2026

*or while stocks last*

Ref. Code: TOMY-2025

1800 008 453 | [sales.au@bio-strategy.dksh.com](mailto:sales.au@bio-strategy.dksh.com) | [www.bio-strategy.com](http://www.bio-strategy.com) | [labshop.dksh.com.au](http://labshop.dksh.com.au)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Delivering Growth

