


# Coverslips for Removable Chambers

**In Stock**

**Only \$89**

Usually \$183.50

Code: 10811



Coverslips for immunofluorescence assays using the ibidi Removable Chambers

- High quality #1.5H glass for mounting and high-end imaging of immunofluorescence stainings
- 24mm x 60mm size for perfect fit with the ibidi Removable Chambers
- D 263 M Schott Glass Material
- No. 1.5H, 170µm (+/-5µm) Thickness
- No autofluorescence and Low birefringence
- Non-Sterile

**Ideal for the 3 Well, 8 Well, and 12 Well Chamber, removable**



Valid Until: 31 Mar 2026

Ref. Code: IBID-10811

0800 34 24 66 | [sales.nz@bio-strategy.dksh.com](mailto:sales.nz@bio-strategy.dksh.com) | [www.bio-strategy.com](http://www.bio-strategy.com) | [labshop.dksh.co.nz](http://labshop.dksh.co.nz)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Delivering Growth