

# Epredia Slides and Coverslips

POWERED BY SHANDON  MICROM  MENZEL-GLÄSER  Richard-Allan Scientific 

## SuperFrost Slides


- SuperFrost slides have a coloured writing area on one side
- A range of colours allows specimens to be labelled and archived simply and clearly
- The writing area is resistant to all common solvents



### Cut Edge Slides

- Economical solution for routine applications and manual work
- Not suitable for use in automatic Printers

Description	Pack Size	Code
Slide SuperFrost Cut Edge White 50pcs	40 Pack	EPBRSF21201 ●
Slide SuperFrost Cut Edge Blue 50pcs	40 Pack	EPBRSF21202 ●
Slide SuperFrost Cut Edge Pink 50pcs	40 Pack	EPBRSF21203
Slide SuperFrost Cut Edge Yellow 50pcs	40 Pack	EPBRSF21204
Slide SuperFrost Cut Edge Green 50pcs	40 Pack	EPBRSF21205
Slide SuperFrost Cut Edge Orange 50pcs	40 Pack	EPBRSF21206



● = Top Seller

### Ground Edge 90° Slides

- Ground edges are polished and help prevent accidental cuts
- Suitable for automatic Printers
- The large surface simplifies picking up of the slides with gripper

Description	Pack Size	Code
Slide SuperFrost Grd Edge 90 Deg Green 50pcs	40 Pack	EPBRSF41290 ●
Slide SuperFrost Grd Edge 90 Deg Orange 50pcs	40 Pack	EPBRSF41291
Slide SuperFrost Grd Edge 90 Deg Purple 50pcs	40 Pack	EPBRSF41295 ●
Slide SuperFrost Grd Edge 90 Deg White 50pcs	40 Pack	EPBRSF41296 ●
Slide SuperFrost Grd Edge 90 Deg Blue 50pcs	40 Pack	EPBRSF41297 ●
Slide SuperFrost Grd Edge 90 Deg Pink 50pcs	40 Pack	EPBRSF41298 ●
Slide SuperFrost Grd Edge 90 Deg Yellow 50pcs	40 Pack	EPBRSF41299 ●



● = Top Seller

### Ground Edge 45° Slides

- Ground edges are polished and help prevent accidental cuts
- Highly suitable for automatic Printers
- Ensures uniform blood smearing as the slide can be placed immediately at the correct angle

Description	Pack Size	Code
Slide SuperFrost Grd Edge 45 Deg White 50pcs	40 Pack	EPBRSF41206 ●
Slide SuperFrost Grd Edge 45 Deg Blue 50pcs	40 Pack	EPBRSF41207
Slide SuperFrost Grd Edge 45 Deg Pink 50pcs	40 Pack	EPBRSF41208
Slide SuperFrost Grd Edge 45 Deg Yellow 50pcs	40 Pack	EPBRSF41209
Slide SuperFrost Grd Edge 45 Deg Green 50pcs	40 Pack	EPBRSF41210



● = Top Seller

POWERED BY

SHANDON

MICROM

MENZEL-GLÄSER

Richard-Allan  
Scientific

## Plain and Twin Frost Slides

- Made from extra-white soda-lime glass
- Clean, washed and polished for immediate use
- Vacuum packing protects the glass against moisture and climate-related damage to slides

### Plain Slides

Description	Pack Size	Code
<b>Plain Slides</b>		
Slide Plain Grd Edge 90 Deg 50pcs	40 Pack	EPBRS31013A
Slide Plain Grd Edge 45 Deg 50pcs	40 Pack	EPBRS31103A

### Twin Frost Slides

Description	Pack Size	Code
Slide Twin Frost Cut Edge 50pcs	40 Pack	EPBRS21102A40 ●
Slide Twin Frost Grd Edge 90 Deg 50pcs	40 Pack	EPBRS41014A ●
Slide Twin Frost Grd Edge 45 Deg 50pcs	40 Pack	EPBRS41104A/40

● = Top Seller

## Adhesion Slides

- Reliable and uniform coating
- Minimal sample loss

### SuperFrost Plus Slides

Electrostatically charged adhesion slides with excellent adhesive properties for high reliability in instrument-based applications.

- For tissue slices 2-5µm thick
- Eliminates background staining in standard H&E stains
- No background staining occurs when immunoperoxidase or in-situ DNA procedures are performed

Description	Pack Size	Code
Slide SuperFrost Plus White 72pcs	1 Pack	EPBRSF41296SP ●

*Other colours are available upon request*

● = Top Seller

### SuperFrost Polysine Slides




- Polysine offers adhesion that is superior to glue, protein or silane treated slides
- For tissue slices 2-5µm thick
- Electrostatically and biochemically adhesive slides
- The preparation is first adhered by electrostatic attraction and then fixed in place by biochemical binding.
- Polysine adhesion slides are best suited to paraffin-embedded tissue slices fixed in formalin, alcohol or Bouins solution

Description	Pack Size	Code
Slide SuperFrost Polysine Coated White 72pcs	1 Pack	EPBRSF41296PL

### SuperFrost Plus Gold Slides

- Revolutionary glass adhesion technology that first attracts and then chemically bonds fresh or frozen tissue slices to their surface
- Particularly suited to poorly adhesive tissue samples (such as bone, breast and brain) that are larger than 5µm
- During immunohistochemical analysis of estrogen and progesterone receptors in fresh frozen breast tissue, SuperFrost Plus Gold slides assist the tissue adherence

Description	Pack Size	Code
Slide SuperFrost Plus Gold 72pcs	1 Pack	EPBRSF41296SPG



POWERED BY SHANDON MICROM MENZEL-GLÄSER Richard-Allan Scientific

## Coverslips

- Archiving for eternity
- Highly transparent and achromatic – practically free of any colour distortion
- High optical homogeneity – practically free of inclusions, streaks and bubbles
- High spectral transmission; Good planarity; Very good chemical stability
- Manufactured from D 263 M borosilicate glass from Schott AG

### Coverslips Standard

Description	Pack Size	Code
Coverslip 18x18mm No. 1 100pcs	10 Pack	EPBRCS18181GP
Coverslip 22x22mm No. 1 100pcs	10 Pack	EPBRCS22221GP ●
Coverslip 22x32mm No. 1 100pcs	10 Pack	EPBRCS22321GP
Coverslip 22x40mm No. 1 100pcs	10 Pack	EPBRCS22401GP ●
Coverslip 22x50mm No. 1 100pcs	10 Pack	EPBRCS22501GP
Coverslip 22x60mm No. 1 100pcs	10 Pack	EPBRCS22601GP
Coverslip 22x70mm No. 1 100pcs	10 Pack	EPBRCS22701GP
Coverslip 24x24mm No. 1 100pcs	10 Pack	EPBRCS24241GP
Coverslip 24x32mm No. 1 100pcs	10 Pack	EPBRCS24321GP
Coverslip 24x40mm No. 1 100pcs	10 Pack	EPBRCS24401GP
Coverslip 24x50mm No. 1 100pcs	10 Pack	EPBRCS24501GP ●
Coverslip 24x60mm No. 1 100pcs	10 Pack	EPBRCS24601GP
Coverslip 50x76mm No. 1 100pcs	10 Pack	EPBRCS50761GP

● = Top Seller

### Coverslips Circle

Description	Pack Size	Code
Coverslip Circle 10mm No.1 100pcs	10 Pack	EPBRCS101GP
Coverslip Circle 12mm No.1 100pcs	10 Pack	EPBRCS121GP
Coverslip Circle 13mm No.1 100pcs	10 Pack	EPBRCS131GP
Coverslip Circle 18mm No.1 100pcs	10 Pack	EPBRCS181GP
Coverslip Circle 25mm No.1 100pcs	10 Pack	EPBRCS251GP
Coverslip Circle 28mm No.1 100pcs	10 Pack	EPBRCS281GP

### Coverslips for Automated Coverslippers

Description	Pack Size	Code
Coverslip 24x50mm No. 1 100pcs	10 Pack	EPBRCSV24501GPA

Other sizes and thicknesses (0,1,1.5 and 2) are available upon request.

Thickness	
No. 0	0.08-0.12mm
No. 1	0.13-0.16mm
No. 1.5	0.16-0.19mm
No. 2	0.19-0.23mm



Contact us for more information  
and to place your orders

**Bio-Strategy Pty Ltd.**

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.