

TOMY New Long Basket!



Save You Lots

Easy loading/
unloading

Excellent draining
efficiency

Smoothly
unloading articles

Punched-out Design

for conducting steam well.

**1. Lightweight Design
with Easy-to-hold Handle**

2. No Mesh Part
for draining without spilling liquid

3. Small Punch Hole Design
prevent articles from getting
stuck in holes.

4. No Mesh in the Bottom
no liquid leakage even if sterilizing
bag is broken.




Goes Directly into the Chamber.

Easy loading/unloading.



As a Waste Bin

Can be used as a waste
bin with sterilizing bag.
The lid to reduce the
smell is also available.



* Do not use the lid for fine
steam circulation.

Long Baskets

	Applicable Autoclaves		Model	Inside dimension (diameter×height)	Applicable Sterilizing Bags
	SX-300 • SX-300E	Basket	ACA-310TP	φ 300×311mm	DAB-321, DAB-60L
		Lid	300TL	—	
	SX-500 • SX-500E	Basket	ACA-325TP	φ 300×491mm	DAB-321, DAB-60L
		Lid	300TL	—	
	SX-700 • SX-700E	Basket	ACA-710TP	φ 345×527mm	DAB-700, DAB-60L
		Lid	700TL	—	
	FLS-1000	Basket	ACA-1000TP	φ 415×456mm	—
		Lid	1000TL	—	

Sterilizing Bags

	Applicable Autoclaves	Model	Max Allowable Temperature	Dimension(mm)	Remarks
	ES-215	DAB-40M	132℃	405W×660Hmm, 45μm	Polypropylene, 100pcs/set
	ES-315 SX-300 • SX-300E SX-500 • SX-500E	DAB-60L	132℃	610W×810Hmm, 45μm	Polypropylene, 100pcs/set
		DAB-321	140℃	550W×750Hmm, 40μm	Polyester, 100pcs/set
	SX-700 • SX-700E	DAB-60L	132℃	610W×810Hmm, 45μm	Polypropylene, 100pcs/set
		DAB-700	140℃	600W×800Hmm, 40μm	Polyester, 100pcs/set


■ Baskets, Buckets

	Autoclave	Inside dimension (diameter × height)	Model	Remarks
Stainless Basket 	ES-215	φ224 × 376mm	ACA-215B	1pc installable
	ES-215/SR-240	φ224 × 181mm	ACA-240B	2pcs installable
	ES-315 SX-300 • SX-300E SX-500 • SX-500E SDL-320	φ300 × 181mm	ACA-315B	3pcs installable (2pcs installable for SX-300/SX-300E)
	SX-700 • SX-700E	φ345 × 181mm	ACA-700B	3pcs installable
	FLS-1000	φ415 × 233mm	ACA-1000B	2pcs installable
Stainless Bucket 	ES-215	φ222 × 364mm	ACA-215K	1pc installable
	SX-300 • SX-300E SX-500 • SX-500E	φ299 × 279mm	ACA-300K	Without hook for stacking
	ES-315 SX-300 • SX-300E SX-500 • SX-500E SDL-320	φ299 × 279mm	ACA-315K	With hook for stacking, 2pcs installable (1pc installable for SX-300/SX-300E)
	ES-315 SX-500 • SX-500E	φ299 × 519mm	ACA-500K	
	SX-700 • SX-700E	φ344 × 279mm	ACA-700K	With hook for stacking, 2pcs installable
			ACA-701K	Without hook for stacking
	FLS-1000	φ415 × 239mm	ACA-1000K	1pc installable

■ Sterilizing Container for Wastes

Model	Specifications	Remarks
ST-ZERO 	Material: Stainless Steel(SUS304) Dimensions(mm): 106W × 350D × 80H Weight: 1.2kg Accommodate up to 25 x 25 ml disposable pipettes	Two stands included

Designed non-preventing air flow in Biosafety cabinet



Leakage prevention design for vertical loading



Sales Office:  Digital Biology®

TOMY DIGITAL BIOLOGY CO., LTD.

3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan

e-mail: info@digital-biology.co.jp

URL: <http://www.digital-biology.co.jp>

phone: +81-3-5971-8160 fax: +81-3-3970-6036

TOMY SEIKO CO., LTD.