

Pall Multi-well Plate Vacuum Manifold



Ref. Code:
PALL-FOLD

Free Demo!

Quick, simple and efficient

The multi-well plate vacuum manifold is an anodized aluminium manifold that has been designed and optimised for the vacuum filtration of AcroPrep, AcroPrep Advance, and AcroWell multi-well filter plates.

- Comes complete with the necessary O-ring and gasket. The control block includes the vacuum pressure gauge, vacuum metering valve, vacuum release valve, and 1/4 in. hose barb for vacuum line attachment
- Vacuum manifold unit includes a Delrins spacer block designed to accommodate standard 350µL receiver plates. The spacer block has been optimised to reduce the space between the receiver plate and the filter plate during vacuum filtration
- Optional spacer block available for use with 1 ml receiver plates
- Adapter collar holds filter plates tightly to receiver plates for centrifugation

Description	Code	List Price (excl. GST)
Multi-well Plate Manifold	PALL5017	\$1,574.00



**FREE
DEMO!**

BIO-STRATEGY PTY LTD.

T: 1800 008 453
W: bio-strategy.com
E: sales.au@bio-strategy.com



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.