

Corning® Falcon® Tubes and Strainers

Simplify your Cell Workflow

Corning® Falcon® Conical Centrifuge Tubes are easy-to-use, high-quality tools to protect your precious samples. State-of-the-art design and manufacturing create tubes engineered to provide high strength, wide temperature stability, and to perform in your critical applications. Corning® Falcon® Conical Centrifuge Tubes are ideal for cell centrifugation, pelleting and separation by density gradients.

- Sterilised by gamma irradiation, noncytotoxic and nonpyrogenic
- Dark blue printed graduations and white writing patch
- Polyethylene flat-top screw cap allows one hand manipulation and provides a level writing area

Falcon® cell strainers are a sterile, easy-to-use device for rapidly isolating primary cells to consistently obtain a uniform single-cell suspension from tissues. Applications include stem cells and primary cells.

- A faster and easier alternative to gauze filtration in procedures involving dissociation of cells from either clumps or primary tissues
- Consistently obtain a more uniform single-cell suspension
- Three nylon mesh pore sizes for optimal performance in a variety of applications
- Sterilised and conveniently accessible in individual packaging
- Extended lip on strainer enables aseptic handling with forceps
- Strainers are made of a strong nylon mesh with 40, 70, or 100µm pores that are evenly spaced for consistent results
- Moulded color-coded polypropylene frame with tab enables easy handling
- Sterile
- Individually packaged

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

VALID UNTIL:
31 March 2019

Ref. Code:
BDAA-FLOW

**Save
Now!**



Description	Pack Size	Availability	Code	Sale Price (excl. GST)
Falcon® 50ml High Clarity PP Centrifuge Tube, Conical Bottom, Sterile	500 items	stocked	BDAA352070	\$99.00
Falcon® 15ml High Clarity PP Centrifuge Tube, Conical Bottom, with Dome Seal Screw Cap, Sterile	500 items	stocked	BDAA352096	\$93.00
Falcon® 40µm Cell Strainer, Blue, Sterile, Individually Packaged	50 items	stocked	BDAA352340	\$135.00
Falcon® 70µm Cell Strainer, White, Sterile, Individually Packaged	50 items	stocked	BDAA352350	\$136.00
Falcon® 100µm Cell Strainer, Yellow, Sterile, Individually Packaged	50 items	stocked	BDAA352360	\$140.00

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.