

# 15% Off All Heidolph Products

## Heidolph Hei-Vap Intelligent Evaporation

### Ask for a demo today

Whether you wish to distil small samples for analysis, prepare large batches, or continuous evaporation – the Hei-VAP range has the right solution for every task. For maximum efficiency and precision, we recommend using our complete solutions, including peripheral devices, which are configurable to your application.

**Leading Safety Standards** – All models come with a safety heating bath with handles, pour spout, and cable coupling with protection class IP 67 to prevent short-circuiting and corrosion. The residual heat indicator also alerts you to a temperature above 50°C when in stand-by mode. In the event of continued temperatures exceeding by >5°C, the heating is turned off to prevent running dry.

**Superior Ease of Use** – Simplify your day-to-day work with smart details such as a detachable panel, control knobs for direct control, the Easy-Clip for removing flasks, and the patented vapor tube. All glass accessories such as distillation spiders, evaporating cylinders, or evaporating flasks up to 5L can be used with the universal 5L heating bath without retrofitting.

**Reduced Cost of Ownership** – All components complement each other and are made for years of continuous operation: PTFE vacuum seals, UV resistant tubes with best chemical stability, condensate-repellent varnish, sealed housing, maintenance-free and sparkless motors. FDA compliant materials increase the range of applications without additional accessories.

**Made in Germany** – All Heidolph products are developed and manufactured at our Schwabach headquarters in Nuremberg, where they undergo multi-stage quality checks in development and production. Even in continuous operation, our powerful, no-maintenance motors ensure consistent results and prevent downtimes and expensive repairs.

**VALID UNTIL:**  
31 March 2019  
**Ref. Code:**  
HEID-VAP

**15%  
OFF**



**Hei-VAP series**  
Up to 5 litre flask distillation

**Hei-VAP Industrial series**  
Up to 20 litre flask distillation

**Hei-PLATE series**  
Heating and stirring

**Hei-MIX series**  
Vortexing and mixing

**Hei-FLOW series**  
Dosing and dispensing

**Hei-TORQUE series**  
Overhead stirring

Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.