

Miele PG8583 Laboratory Washer Package 1

VALID UNTIL:
31 March 2019
Ref. Code:
MIEL-G8583

**10% Package Discount and free commissioning!
(city metro areas only)**

Miele Professional PG8583 have been the standard in Laboratory Underbench/ Free-Standing Washers for the past three years. Their stainless-steel interior combining powerful pump with integrated heating elements enable safe washing cycles up to 93°C. The EcoDry passive drying function automatically opens the door at the end of the run when temperature has dropped below 70°C allowing the load to dry and cool down faster. The unit also features spray arm and spray pressure monitoring with alerts to advise you if a blockage or pressure deviation has occurred.

- 15 programs including new "Mini" program for lightly soiled glassware
- ADP model including pump for AD water and waterproofing kit
- Includes upper and lower baskets and inserts/modules for washing wide neck and narrow neck glassware
- Dispensing pumps for liquid detergent and liquid neutraliser
- Compact: 600x600x850mm (820mm with lid removed)
- Many other options and accessories to select from

Versatile washing solution from a German manufacturer – Miele

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com



**ON
SALE FOR
\$11,653.00**

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
Washer PG8583 + Accessories for washing wide neck and narrow neck glassware at the same time	MIELPG8583-Package-1	\$12,948.00	\$11,653.00
Commissioning/Training Australian Metro Areas	MIELMISC		No Charge

Please note: Price includes washer, waterproofing kit, A102 Upper Basket with 2xAK12 inserts for wide neck glassware, A150 Lower Basket for modules and 2xA301/2 modules (each with 3x6 = 18 injector nozzles; 160mm high), customer to organise local plumber & electrician for installation and cover local freight. Commissioning to be done by Miele Professional once unit is installed (included in price, city metro area's only).

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.