

Pall Centrifugal Devices

Save on Spin Devices

These devices are used for simple, reliable concentrating and desalting of samples.

Select the low protein binding Omega™ membrane for ultrafiltration applications. The Bio-Inert® membrane is ideal for microfiltration of biological solutions. The GHP membrane is ideal for sample preparation in analytical chemistry (HPLC, GC, IC and MS). It is suitable for aqueous, organic and aggressive organic samples and for rapid processing of samples.

- Rapid processing of samples
- Typical recoveries are greater than 90%
- Housing: Low binding PP
- Ultrasonically welded seals prevent bypass or seal failure

**10%
OFF**

**VALID UNTIL:
31 March 2019**

**Ref. Code:
PALL-CENT**



Description	Pack Size	Code	Sale Price (excl. GST)
Nanosep® MF with Bio-Inert® membrane, 0,45 µm, Purple	24 items	PALLODM45C33	\$78.30
Nanosep® MF with GHP membrane, 0,45 µm, Clear	500 items	PALLODGHPC35	\$1,420.20
Nanosep® with Omega™ membrane, 30 kD, Red	24 items	PALLOD030C33	\$106.20
Nanosep® with Omega™ membrane, 100 kD, Clear	24 items	PALLOD100C33	\$107.10
Nanosep® MF with GHP membrane, 0,45 µm, Clear	100 items	PALLODGHPC34	\$378.00
Nanosep® MF with Bio-Inert® membrane, 0,2 µm, Turquoise	100 items	PALLODM02C34	\$273.60
Nanosep® MF with Bio-Inert® membrane, 0,45 µm, Purple	100 items	PALLODM45C34	\$273.60
Nanosep® with Omega™ membrane, 10 kD, Blue	24 items	PALLOD010C33	\$110.70
Nanosep® MF with Bio-Inert® membrane, 0,2 µm, Turquoise	500 items	PALLODM02C35	\$1,359.90
Nanosep® with Omega™ membrane, 300 kD, Orange	100 items	PALLOD300C34	\$357.30
Nanosep® with Omega™ membrane, 30 kD, Red	500 items	PALLOD030C35	\$1,494.90
Nanosep® with Omega™ membrane, 300 kD, Orange	500 items	PALLOD300C35	\$1,494.90
Nanosep® with Omega™ membrane, 10 kD, Blue	100 items	PALLOD010C34	\$478.80
Nanosep® with Omega™ membrane, 10 kD, Blue	500 items	PALLOD010C35	\$1,497.60
Nanosep® with Omega™ membrane, 300 kD, Orange	24 items	PALLOD300C33	\$107.10
Nanosep® MF with Bio-Inert® membrane, 0,2 µm, Turquoise	24 items	PALLODM02C33	\$78.30
Nanosep® MF with Bio-Inert® membrane, 0,45 µm, Purple	500 items	PALLODM45C35	\$1,170.00
Nanosep® with Omega™ membrane, 3 kD, Grey	500 items	PALLOD003C35	\$1,494.90
Nanosep® with Omega™ membrane, 100 kD, Clear	500 items	PALLOD100C35	\$1,494.90
Nanosep® with Omega™ membrane, 3 kD, Grey	100 items	PALLOD003C34	\$357.30
Nanosep® with Omega™ membrane, 100 kD, Clear	100 items	PALLOD100C34	\$357.30
Nanosep® with Omega™ membrane, 3 kD, Grey	24 items	PALLOD003C33	\$106.20

More codes on the next page...

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

Pall Centrifugal Devices

Description	Pack Size	Code	Sale Price (excl. GST)
Microsep® with Omega™ Membrane, 3 kD, Grey	24 items	PALLMCP003C41	\$139.83
Microsep® with Omega™ Membrane, 10 kD, Blue	24 items	PALLMCP010C41	\$139.83
Microsep® with Omega™ Membrane, 30 kD, Red	24 items	PALLMCP030C41	\$139.83
Microsep® with Omega™ Membrane, 100 kD, Clear	24 items	PALLMCP100C41	\$139.83
Microsep® with Omega™ Membrane, 0.2 µm, Aqua	24 items	PALLMCPM02C67	\$139.83
Microsep® with Omega™ Membrane, 0.45 µm, Wildberry	24 items	PALLMCPM45C67	\$139.83
Microsep® with Omega™ Membrane, 10 kD, Blue	100 items	PALLMCP010C46	\$502.68
Microsep® with Omega™ Membrane, 30 kD, Red	100 items	PALLMCP030C46	\$502.68
Microsep® with Omega™ Membrane, 100 kD, Clear	100 items	PALLMCP100C46	\$502.68
Microsep® with Omega™ Membrane, 1 kD, Yellow	100 items	PALLMCP001C46	\$507.99
Microsep® with Omega™ Membrane, 1 kD, Yellow	24 items	PALLMCP001C41	\$141.60
Microsep® with Omega™ Membrane, 0.2 µm, Aqua	100 items	PALLMCPM02C68	\$507.99
Microsep® with Omega™ Membrane, 3 kD, Gray	100 items	PALLMCP003C46	\$502.68
Microsep® with Omega™ Membrane, 0.45 µm, Wildberry	100 items	PALLMCPM45C68	\$502.68
Macrosep® with Omega™ Membrane, 3 kD, Gray	6 items	PALLMAP003C36	\$73.46
Macrosep® with Omega™ Membrane, 10 kD, Blue	6 items	PALLMAP010C36	\$73.46
Macrosep® with Omega™ Membrane, 30 kD, Red	6 items	PALLMAP030C36	\$73.46
Macrosep® with Omega™ Membrane, 100 kD, Clear	6 items	PALLMAP100C36	\$76.11
Macrosep® with Omega™ Membrane, 3 kD, Gray	24 items	PALLMAP003C37	\$260.19
Macrosep® with Omega™ Membrane, 10 kD, Blue	24 items	PALLMAP010C37	\$260.19
Macrosep® with Omega™ Membrane, 100 kD, Clear	24 items	PALLMAP100C37	\$260.19
Macrosep® with Omega™ Membrane, 3 kD, Gray	100 items	PALLMAP003C38	\$1,015.98
Macrosep® with Omega™ Membrane, 10 kD, Blue	100 items	PALLMAP010C38	\$1,015.98
Macrosep® with Omega™ Membrane, 30 kD, Red	100 items	PALLMAP030C38	\$1,015.98
Macrosep® with Omega™ Membrane, 0.2 µm, Aqua	100 items	PALLMAPM02C68	\$1,015.98
Macrosep® with Omega™ Membrane, 0.45 µm, Wildberry	100 items	PALLMAPM45C68	\$1,015.98
Macrosep® with Omega™ Membrane, 1 kD, Yellow	6 items	PALLMAP001C36	\$74.34
Macrosep® with Omega™ Membrane, 30 kD, Red	24 items	PALLMAP030C37	\$260.19
Macrosep® with Omega™ Membrane, 0.2 µm, Aqua	24 items	PALLMAPM02C67	\$260.19
Macrosep® with Omega™ Membrane, 1 kD, Yellow	100 items	PALLMAP001C38	\$1,065.54
Macrosep® with Omega™ Membrane, 1 kD, Yellow	24 items	PALLMAP001C37	\$262.85
Macrosep® with Omega™ Membrane, 100 kD, Clear	100 items	PALLMAP100C38	\$1,015.98

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

<https://shop.bio-strategy.com>



bio-strategy
delivering technology

Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.