

PHCbi TwinGuard® -86°C Ultra-Low Temperature (ULT) Freezer

Uncompromising storage for your samples!

The new PHCbi TwinGuard® -86°C ultra-low temperature (ULT) freezer satisfies an industry demand for safe, long-term storage for highly-valued samples. Two independent refrigeration systems combined with optional liquid N2 or liquid CO2 back-up systems offer levels of protection that is unmatched in the marketplace.

- Unique independent dual cooling system - maintain -65°C with a single compressor
- Improved 5-fold Gasket – increases insulation and reduces frosting within the chamber
- USB accessible on board datalogging and diagnostics
- Filterless design – no filters to clean and change
- VIP+ Insulation – store more samples than ever before within the same footprint
- Various additions available – CO2 backup, flexible inner doors, IQ/OQ, racks and cryoboxes, external monitoring, service and maintenance options

PHCbi TwinGuard is available in 3 upright sizes. Units in Stock. Contact your Bio-Strategy sales representative for more information and TwinGuard your samples today!

Description	Code
PHCbi TwinGuard ULT, 360L, WxDxH 670x882x1840mm	SANYMDF-DU302VX-PE
PHCbi TwinGuard ULT, 528L, WxDxH 790x882x1993mm	SANYMDF-DU502VX-PE
PHCbi TwinGuard ULT, 729L, WxDxH 1030x882x1993mm	SANYMDF-DU702VX-PE
Accessories	
Liquid CO2 Backup	SANYMDF-UB7-PW
Set of 5 inner doors, includes 1 shelf for DU502VX	SANYMDF-5ID5-PW
Set of 4 inner doors for DU502VX	SANYMDF-5ID4-PW
Set of 5 inner doors for DU702VX, includes additional shelf	SANYMDF-7ID5-PW
Set of 4 inner doors for DU702VX	SANYMDF-7ID4-PW

Safeguard your samples with TwinGuard!



TwinGuard

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.