

LabAdvantage Barrier Pipette Tips only \$52.50

Take advantage of this special price offer!

All products are manufactured from the finest raw materials using precise production techniques that provide consistent quality performance. All primary products are made from virgin polymers and contain no recycled, reprocessed, or altered plastics that might introduce contamination or defects.

Every lot of each product can be traced to the corresponding raw materials and production date. In addition to continual in-process control monitoring, all product lots are sampled after packaging and tested in our quality control laboratory as well as independent laboratories. The result is a contaminant-free product of the highest quality.

- Ideal for molecular biology applications like PCR
- Features an inert, hydrophobic barrier that prevents pipettor contamination
- Aerosol barrier will not trap liquid or contaminate samples in the event of over-pipetting
- Universal fit for research-grade pipettors
- Available with Low Binding Surface Technology

VALID UNTIL:
31 March 2019
Ref. Code:
LABAD-TIPS

**Only
\$52.50!**

Description	Unit of Measure	Code	Sale Price (excl. GST)
10µL "E" Barrier Tip, Sterile, Low Binding, Extended Length, 10x96 Tips (Racked)	960 items	SORELAC15024T	\$52.50
10µL Barrier Tip, Sterile, Low Binding, 10x96 Tips (Racked)	960 items	SORELAC28204T	\$52.50
10µL Barrier Tip, Sterile, Low Binding, Extended Length, 10x96 Tips (Racked)	960 items	SORELAC30364T	\$52.50
20µL Barrier Tip, Sterile, Low Binding, 10x96 Tips (Racked)	960 items	SORELAC30374T	\$52.50
200µL Barrier Tip, Sterile, Low Binding, 10x96 Tips (Racked)	960 items	SORELAC30554T	\$52.50
200µL Barrier Tip, Sterile, Low Binding, Multi-fit, 10x96 Tips (Racked)	960 items	SORELAC17374T	\$52.50
1000µL Barrier Tip, Sterile, Low Binding, 10x100 Tips (Racked)	1000 items	SORELAC14204T	\$52.50
1000µL Barrier Tip, Sterile, Low Binding, Extended Length, 8x96 Tips (Racked)	768 items	SORELAC34004	\$52.50



Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.