

Tomy SX Series Top Loading Autoclave Promotion

VALID UNTIL:
30 March 2019

Ref. Code:
TOMY-SX

Including 1 Free Preventative Maintenance Inspection!

Convenient sterilising in a compact installation space
Three Capacities – choose the one that's right for you

- Three operating temperature ranges: warming 45 to 95°C, heating 45 to 104°C and sterilising 105 to 135°C
- Safety devices: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection and safety valve
- Integrated (but removable) condensate bottle to trap condensate from chamber during venting
- Simple lid opening using either hand or foot
- Cooling fan to reduce run time
- 1 Free Preventative Maintenance visit (metro area's only)
- Supplied with 1 basket (SX-300E) or 2 baskets (SX-500E & SX-700E)
- Digital temperature and time display with analogue gauge for pressure
- Five sterilisation profiles to select from (in addition to normal sterilisation)

Speak to your local Bio-Strategy Equipment Specialist for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

ON SALE

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
Top Loading 44 Litre Autoclave	TOMYTMYSX-300E	\$11,548.00	\$8,400.00
Top Loading 58 Litre Autoclave	TOMYTMYSX-500E	\$11,654.00	\$8,500.00*
Top Loading 79 Litre Autoclave	TOMYTMYSX-700E	\$13,725.00	\$9,975.00*

**Ex stock*



Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.