

The Most Economic VWR Quality Pipettors EXCLUSIVE TO BIO-STRATEGY

Time to replace your old pipettors!

We offer 3-lines of VWR pipettors at special prices this Summer. All pipettors are autoclavable, compatible to most pipette tips, UV resistant and come with a 3-year warranty.

Stock availability can vary. Contact our Customer Service on 1800 008 453 or email sales.au@bio-strategy.com for details.

VWR Ultra-High-Performance (UHP) Pipettors

VWR Ultra-High-Performance pipettes offer superior accuracy, precision and ergonomic design. Lightweight and durable. Fully autoclavable. Single and multichannel models. 3-year warranty. The VWR premium pipettor!

- Very low plunger forces and near universal tip capability
- Automatic locking volume adjustment wheel and levered colour coded ejection

Description	Accuracy	Imprecision	Availability	Code	Sale Price (excl. GST)
Pipettor Ultra-High-Performance					
0.2-2.0µL	±12.0 - ±1.5 %	<6.0 - <0.7 %	enquire	VWRI613-1488	\$369.00
2-20µL	±3.0 - ±0.8 %	<1.5 - <0.3 %	enquire	VWRI613-1490	\$375.00
20-200µL	±1.2 - ±0.6 %	<0.8 - <0.2 %	stocked	VWRI613-1492	\$369.00
100-1000µL	±0.9 - ±0.6 %	<0.6 - <0.2 %	stocked	VWRI613-1493	\$369.00
1.0-10µL	±8.0 - ±2.0 %	<6.0 - <1.2 %	stocked	VWRI613-1494	\$799.00
5-50µL	±4.0 - ±1.6 %	<2.5 - <0.6 %	stocked	VWRI613-1495	\$799.00
8-Channel, 20-200µL	±3.0 - ±1.0 %	<1.5 - <0.6 %	enquire	VWRI613-1496	\$864.00
12-Channel, 1.0-10µL	±8.0 - ±2.0 %	<6.0 - <1.2 %	enquire	VWRI613-1498	\$1,143.00
12-Channel, 20-200µL	±3.0 - ±1.0 %	<1.5 - <0.6 %	stocked	VWRI613-1500	\$989.00
VWR UHP Packs					
VWR UHP 5 Pack incl. 5 pipettes (0.1-2, 0.5-10, 2-20, 20-200 and 100-1000µl), 3 racks of tips (10, 200 and 1000µl) and 8 pipette multiple stand			enquire	VWRI613-5464	\$1,791.00
VWR UHP 3 Pack incl. 3 pipettes (2-20, 20-200, 100-1000µl), 2 racks of tips (200 and 1000µl) and 8 pipette multiple stand			enquire	VWRI613-5467	\$1,057.00

VWR Standard Line Economic Pipettors

These lightweight Standard Line pipettes provide accuracy and precision. Ergonomic features include contoured handle shape and a soft spring system. 3-year warranty, autoclavable.

- Easily adjustable ejector height system for many tips
- UV resistant and easy to calibrate
- Colour coded shaft for easy volume identification

Description	Accuracy	Imprecision	Availability	Code	Sale Price (excl. GST)
Pipettor Standard Line					
0.2-2µL	±12.0 - ±1.5 %	≤6.0 - ≤0.7 %	stocked	VWRI613-2891	\$199.00
10-100µL	±1.6 - ±0.8 %	≤0.8 - ≤0.2 %	enquire	VWRI613-2894	\$199.00
20-200µL	±1.2 - ±0.6 %	≤0.6 - ≤0.2 %	enquire	VWRI613-2895	\$199.00
100-1000µL	±1.6 - ±0.6 %	≤0.4 - ≤0.15 %	enquire	VWRI613-2896	\$199.00
500-5000µL	±1.2 - ±0.5 %	≤0.5 - ≤0.15 %	enquire	VWRI613-2897	\$199.00
1000-10000µL	±2.5 - ±0.5 %	≤0.6 - ≤0.2 %	enquire	VWRI613-2898	\$199.00
Pipettor Standard Line Pack					
3 Pack (2-20, 20-200, 100-1000µl)			stocked	VWRI613-5469	\$549.00

VALID UNTIL:
31 March 2019
Ref. Code:
VWRI-PIPE

**Special
Prices on
3-lines of
pipettors!**



**ON
SALE**



**ON
SALE**

VWR Ergonomic High Performance (EHP) Pipettors

The Ergonomic High-Performance pipettes provide accuracy and precision you can rely on. 3-year warranty, UV-resistant, autoclavable.

- Significantly reduced plunger forces
- Effortless volume adjustment with lock
- Tapered shaft and adjustable 'Easy Grip' ejector for near universal tip fit

Description	Accuracy	Imprecision	Availability	Code	Sale Price (excl. GST)
Pipettor Ergonomic-High-Performance (EHP)					
0.2-2µL	±12.0 - ±1.5 %	<6.0 - <0.7 %	enquire	VWRI613-5258	\$318.50
0.5-10µL	±1.2 - ±0.6 %	<2.8 - <0.4 %	stocked	VWRI613-5259	\$318.50
2-20µL	±1.0 - ±0.6 %	<1.5 - <0.3 %	enquire	VWRI613-5260	\$318.50
5-50µL	±0.9 - ±0.6 %	<2.0 - <0.4 %	enquire	VWRI613-5261	\$318.50
10-100µL	±0.6 - ±0.5 %	<0.8 - <0.2 %	stocked	VWRI613-5262	\$318.50
20-200µL	±2.5 - ±0.5 %	<0.6 - <0.2 %	stocked	VWRI613-5263	\$318.50
50-250µL	±1.0 - ±0.6 %	<0.4 - <0.3 %	enquire	VWRI613-5264	\$318.50
100-1000µL	±0.9 - ±0.6 %	<0.4 - <0.15 %	stocked	VWRI613-5265	\$318.50
1000-5000µL	±0.6 - ±0.5 %	<0.25 - <0.15 %	stocked	VWRI613-5266	\$324.00
1000-10000µL	±2.5 - ±0.5 %	<0.6 - <0.2 %	enquire	VWRI613-5267	\$324.00
8-Channel, 1-10µL	±8.0 - ±2.0 %	<8.0 - <2.0 %	enquire	VWRI613-5250	\$678.00
8-Channel, 5-50µL	±4.0 - ±1.6 %	<4.0 - <1.6 %	enquire	VWRI613-5251	\$678.00
8-Channel, 20-200µL	±2.0 - ±1.0 %	<1.5 - <0.8 %	enquire	VWRI613-5252	\$678.00
8-Channel, 50-300µL	±1.6 - ±1.0 %	<1.5 - <0.8 %	enquire	VWRI613-5253	\$678.00
Pipettor Ergonomic-High-Performance (EHP) Packs					
3 Pack, incl. 3 pipettes (0.5-10, 10-100 and 100-1000µl), 3 racks of tips (10, 200 and 1000µl) and 4 pipette linear stand			enquire	VWRI613-5466	\$889.00
3 Pack, incl. 3 pipettes (2-20, 20-200, and 100-1000µl), 2 racks of tips (200 and 1000µl) and 4 pipette linear stand			enquire	VWRI613-5465	\$889.00
4 Pack, incl. 4 pipettes (0.5-10, 2-20, 20-200, and 100-1000µl), 3 racks of tips (10, 200 and 1000µl) and 4 pipette linear stand			enquire	VWRI613-5462	\$1,134.00



Special Clearance Sale

Description	Availability	Code	Sale Price (excl. GST)
Pipettor Ergonomic-High-Performance (US)			
8-Channel 1-10µL	While stock lasts	VWRI89079-944	\$789.00
8-Channel 5-50µL	While stock lasts	VWRI89079-946	\$789.00
12-Channel 5-50µL	While stock lasts	VWRI89079-954	\$1,047.00
12-Channel 20-200µL	While stock lasts	VWRI89079-956	\$1,047.00
12-Channel 50-300µL	While stock lasts	VWRI89079-958	\$1,047.00
0.2-2µL	While stock lasts	VWRI89079-960	\$339.00
5-10µL	While stock lasts	VWRI89079-962	\$339.00
5-50µL	While stock lasts	VWRI89079-966	\$317.00
20-200µL	While stock lasts	VWRI89079-970	\$317.00
1-10mL	While stock lasts	VWRI89079-978	\$317.00
Pipette Kit (EHP)			
Pipette Kit 2 - 3 Pipettors (0.2-2, 0.5-10 and 10-100µL)	While stock lasts	VWRI613-5279	\$799.00
Pipette Kit 3 - 3 Pipettors (0.5-10, 10-100 and 100-1000µL)	While stock lasts	VWRI613-5280	\$799.00



Bundle up with VWR Next Generation Pipette Tips for further savings!

Bio-Strategy Pty Ltd.

T: 1800 008 453
 E: sales.au@bio-strategy.com
 www.bio-strategy.com
<https://shop.bio-strategy.com>



bio-strategy
 delivering technology

Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.