

# Ansell's BioClean™ High Quality Solutions for Cleanrooms

**VALID UNTIL:**  
30 September 2019

**Ref. Code:**  
BIOCLEAR

**All BioClean™ (formerly known as Nitritex) products are pristine to ensure products and processes remain protected. All products are processed and packed within an ISO Class 4 or 5 cleanroom.**

Ansell's BioClean™ cleanroom solutions offer a range of consumables including clean and clean/validated sterile disposable gloves, autoclavable and single-use goggles, and disposable garments for both product and personal protection.

#### **Sterile Nitrile Cleanroom Glove products include:**

- BioClean™ Emerald – Wide range of chemical resistance, powder, sulphur and accelerator free, latex-free, textured surface or smooth feel, easy double donning (product code: BENS)
- BioClean™ Indigo - Wide range of chemical resistance, powder-free and latex-free, textured fingers for enhanced grip (product code: BNPLS)
- BioClean™ N-Plus - Wide range of chemical resistance, elbow length protection, powder-free and latex-free, antistatic, non-particulating (product code: BNPS)

#### **Sterile Latex Cleanroom Glove products include:**

- BioClean™ Maxima – Shoulder length protection, powder-free, beaded cuff, flexible and comfortable (product code: BLLS)
- BioClean™ Extra – Elbow length protection, powder-free, beaded cuff, flexible and comfortable (product code: BLAS)

#### **Non-Sterile Latex Cleanroom Glove products include:**

- BioClean™ Legion – Elbow length protection, powder-free, flexible and comfortable (product code: BLA3)

#### **Non-Sterile Nitrile Cleanroom Glove products include:**

- BioClean™ Nerva – Wide range of chemical resistance, elbow length protection, anti-static and beaded cuff (product code: BNAL)
- Omega Nerva™ - Wide range of chemical resistance, elbow length protection, ambidextrous, beaded cuff (product code: ONAW)

#### **Sterile Disposable Garment products include:**

- BioClean-D Coverall with Hood – Elasticated back, hood and ankles, foot-loop to aid smooth zip closure, thumb loops for secure hold, sealable zip cover and antistatic lightweight low-linting Clean Tough material (product code: S-BDCHT)
- BioClean-D Sterile Overboots – Elasticated opening, lightweight low-linting Clean Tough material, secure tie fastenings and slip-resistant sole (product code: S-BDOB)
- BioClean-D Sterile Sleeve Covers – Extra-long length 500mm, antistatic lightweight low-linting Clean Tough material and elasticated openings for firm fit (product code: S-BDSC-L)

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)**

*Product codes and more information on the next page...*



Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

**On Sale,  
Samples  
Available, New  
Stock Arriving  
August**



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

# Ansell's BioClean™ High Quality Solutions for Cleanrooms

## Accessories include:

- BioClean™ Clearview Sterile Visor Face Mask – Anti-fog coating, distortion free visor, latex-free, high bacterial and particle filtration efficiency, looped with blue clip fastener (product code: VFM210-L)
- BioClean™ Clearview Autoclavable Goggles – Super-soft frame with indirect ventilation system, non-linting latex-free head band, toughened polycarbonate, anti-scratch and optically correct lens, tested to withstand >40 autoclave cycles and 12 gamma irradiation cycles (product code: BCAH)
- BioClean™ Clearview Autoclavable Panoramic Goggles – Super-soft frame with indirect ventilation system, ideal for wearing over eye-glasses, non-linting latex-free head band, toughened polycarbonate, anti-scratch and optically correct lens, tested to withstand >40 autoclave cycles (product code: BCAP)
- BioClean™ DB Sterile Pouch-Style Face Mask with Neck Guard – High bacterial and particle filtration efficiency, large breathing chamber, fully enclosed malleable nose-band and extra-long neck guard (product code: DBBS-G)
- BioClean™ Sterile Looped/Tie-On Face Mask – High bacterial and particle filtration efficiency, fully enclosed malleable nose-band, ultrasonically sealed edges and looped with secure fastener or tie-on face mask options available (product code: MEA210-1 / MTA210-1)
- BioClean™ SafeStep NSO Cleanroom Overshoes – Slip resistant, low level particle shedding, durable, elastomer coating and available in 16" (product code: NSO016)



| Description   | Pack Size   | Code         | Sale Price (excl. GST) |
|---|-------------|--------------|------------------------|
| B/Clearview 16CS AC Goggles                         | 12 Items    | NITRBCAH     | <b>\$776.81</b>        |
| B/Clearview 16CSC AC Panoramic Goggles              | 12 Items    | NITRBCAP     | <b>\$769.68</b>        |
| BioClean St.Pouch Mask N.Guard                      | 200 Items   | NITRDBBS-G   | <b>\$567.76</b>        |
| BioClean Emerald St Nit (Sizes Available: 6-10)     | 200 Pairs   | NITRBENS     | <b>\$460.86</b>        |
| BioClean Legion 400 Latex (Sizes Available: S-XL)   | 1,000 Items | NITRBLA301   | <b>\$705.54</b>        |
| BioClean Extra 400 St Latex (Sizes Available: S-XL) | 200 Pairs   | NITRBLAS     | <b>\$612.90</b>        |
| BioClean Maxima St Latex (Sizes Available: 6-10)    | 100 Pairs   | NITRBLLS     | <b>\$1,042.87</b>      |
| BioClean Nerva Nitrile (Sizes Available: S-2XL)     | 1,000 Items | NITRBNAL     | <b>\$741.18</b>        |
| BioClean Indigo St Nit (Sizes Available: 6-10)      | 200 Pairs   | NITRBNPLS    | <b>\$313.57</b>        |
| BioClean N.Plus St Nit (Sizes Available: 6-10)      | 200 Pairs   | NITRBNPS     | <b>\$810.07</b>        |
| BioClean St.Mask E.loop Indv.                       | 600 Items   | NITRMEA210-1 | <b>\$603.39</b>        |
| BioClean St Mask Tie-On Indv.                       | 500 Items   | NITRMTA210-1 | <b>\$620.02</b>        |
| BioClean SafeStep-16 White                          | 400 Items   | NITRNSO016   | <b>\$249.43</b>        |
| Omega Nerva Wh Nit Ambi400 (Sizes Available: S-2XL) | 1,000 Items | NITRONAW     | <b>\$534.50</b>        |
| BioClean-D St.Cvrl Hood (Sizes Available: XS-2XL)   | 20 Items    | NITRS-BDCHT  | <b>\$320.70</b>        |
| BioClean-D St.Cvrl Hood (Sizes Available: 3XL-4XL)  | 15 Items    | NITRS-BDCHT  | <b>\$320.70</b>        |
| BioClean-D St Overboot                              | 75 Pairs    | NITRS-BDOB   | <b>\$382.47</b>        |
| BioClean-D St Sleeve Cvr-500mm                      | 90 Pairs    | NITRS-BDSC-L | <b>\$325.45</b>        |
| BioClean St Loop Visor Facemask                     | 100 Items   | NITRVFM210-L | <b>\$356.33</b>        |

**On Sale  
and Samples  
Available**



Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>



**bio-strategy**  
delivering technology

Prices & discounts above may be subject to change. Not valid with any other special offer or discounts and cannot be applied in retrospect. All offers subject to availability and while stocks last. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.