

The Elga PURELAB Flex 3 for RO and Type I Water

The all-rounder for laboratory-grade water.

Flex 3 system from ELGA Veolia delivers both RO and Ultrapure Type I water from a potable feed source. Its' small and compact design means the unit can be placed on a benchtop or wall mounted.

Features:

- Fitted with 7L reservoir tank
- Real Time TOC monitoring
- Easily maintained with quick and easy consumables to replace
- No nuisance alarms – only warns you when the water quality drops
- Full recirculation through UV and Purification Cartridge – Guaranteed water purity from point of dispense

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

Description	Code
ELGA PURELAB Flex 3 System	ELGAPF3XXXXM1



**VALID UNTIL:
30 September 2019**

**Ref. Code:
ELGA-FLEX3**

**Save
25%!**

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.