

## Grant QBD1 Dry Block Heater

### Robust, compact and with a contemporary design.

The Grant QBD series are a versatile range of high quality digital dry block heating systems with excellent temperature control and a wide choice of interchangeable blocks for maximum flexibility.

The outstanding performance of the QBD series is based on advanced digital temperature control combined with high quality, precision-engineered blocks providing excellent thermal contact. The result is accurate, reproducible and safe heating of your samples.

The QBD1 is a digital system with one optional interchangeable block capacity.

There are an extensive range of optional interchangeable heating blocks to accommodate your sample tubes, from our standard range of blocks or tailored to suit your application.

- Heating range ambient +5°C to 130°C
- Rapid heat-up time
- Stability and uniformity  $\pm 0.1^\circ\text{C}$
- Integral timer function optional external probe for in-sample or in-block temperature control
- Block extraction tool included as standard, for easy and safe removal of blocks
- 2-year warranty

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Digital Dry Block Heater	GRANQBD1	\$1,309.00	\$1,112.00
<b>Best Selling Blocks</b>			
Block 24x1.5ml Tubes	GRANQB-E1	\$205.00	\$175.00
Block 12x16mm $\varnothing$ Tubes	GRANQB-E16	\$205.00	\$150.00

*\*Supplied without optional block.*



**VALID UNTIL:  
30 September 2019**

**Ref. Code:  
GRANQBD1**

**15% off on  
current QBD1  
stock!**

Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.