

## Grant PV-1 Vortex Mixer on Sale

**Exceptionally compact with a low profile and small footprint - fits into almost any location.**

Ideally sized for personal use, the Grant-bio PV-1 vortex mixer is compact, with a low profile and small footprint to sit in almost any location. The PV-1 is so small that it will fit in your hand, yet there is no compromise on performance. The reliable motor produces regulated and reproducible agitation throughout the speed range. It is also extremely quiet.

With variable speed from 750 to 3000rpm, the PV-1 can be used for gentle through to vigorous mixing and pellet resuspension in 0.5ml, 1.5ml and 15ml tubes.

- Adjustable speed control: 750 to 3000rpm
- 'Continuous' or 'touch' operation
- For tubes up to 20mm diameter
- Pressure sensitive cup accommodates tubes up to 20mm diameter
- Reliable and extremely quiet motor produces regulated and reproducible agitation throughout the speed range
- 2-year warranty

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Vortex Mixer	GRANPV-1	\$329.00	<b>\$275.00</b>

*\*Includes 20mm Ø polystyrene mixing cup.*

**Compact  
Affordable  
Vortex  
Mixing!**



Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

**VALID UNTIL:  
30 September 2019**

**Ref. Code:  
GRANPV-1**

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.