

Grant Stirred Waterbath Model T100-ST38

VALID UNTIL:
30 September 2019

Ref. Code:
T100-ST38

38-litre bath with immersion thermostat.

The Grant T100 series of general-purpose heating circulating baths offer accurate temperature control for straightforward applications in the range of 0°C to 100°C. For below ambient temperature control (down to 0°C), an optional cooling device is required. The T100 heating circulator can be used in combination with any Grant stainless steel or plastic tank or with a T-CLAMP for independent use of virtually any vertical sided tank. Our offer has the T100 combined to a 38-litre stainless-steel tank with bridge piece.

- Heating range ambient +5°C to 100°C (below ambient work possible with cooling unit)
- Stability $\pm 0.05^\circ\text{C}$
- Three programmable temperature pre-sets
- Low liquid protection and fixed over temperature cut-out
- Tank external dimensions: 730x335x260mm WxDxH
- 3-year warranty

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Immersion Thermostat	GRANT100	\$1,564.00	\$1,250.00
38-Litre Tank with bridge	GRANST38	\$2,024.00	\$1,620.00
Optional:			
Gabled Hinged SS Lid	GRANSTL38	\$882.00	\$705.00

**Supplied without lid.*

**20%
Package
Discount!**



Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.