

## IKA RCT Basic Stock Sale!

**IKA RCT Basic, the must have hotplate magnetic stirrer for all laboratories.**

- 50-1500rpm, 20L max, 0-310°C
- Integrated Temperature control (incl. PT1000 temperature sensor)
- Digital display of temperature and speed
- Hot-Top indicator – hot surface warning!
- Aluminium heating plate ø135mm, for optimum heat transfer
- Enclosed assembly (IP42)

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)**

**VALID UNTIL:  
30 September 2019**

**Ref. Code:  
IKA-RCT**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
IKA RCT Basic Hotplate Magnetic Stirrer (20003581 – AU Version)	IKAA20002620	<b>\$1,533.00</b>	<b>\$690.00</b>

**Stock  
Sale**



Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.