

# INTEGRA Biosciences Pipetboy Pro Pipette Controller

**Makes liquid handling with serological pipettes accurate, fast and comfortable.**

**Choose your colour: White, Pink, Green, Orange or Blue.**

With over four decades of liquid handling experience and hundreds of thousands of instruments in use worldwide, INTEGRA pipette controllers have become the reference for precise and efficient serological pipetting.

Liquid handling lies at the core of every life science experiment and a key to successful research is the ability to handle liquids with precision and efficiency. The PIPETBOY pro pipetting aid unites all essential characteristics for productive work with serological pipettes: accuracy, speed and ergonomic design.

- Speed range: 2-12ml/s
- LED lighting illuminates your work area
- Pipette compatibility: glass and plastic pipettes from 1 to 100ml
- Filter protection: 0.45µm (included); 0.2µm optional
- Weight: 190g
- 3-year warranty

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
PIPETBOY Pro White	INTG156400	<b>\$819.00</b>	<b>\$475.00</b>
PIPETBOY Pro Pink	INTG156401	<b>\$819.00</b>	<b>\$475.00</b>
PIPETBOY Pro Green	INTG156402	<b>\$819.00</b>	<b>\$475.00</b>
PIPETBOY Pro Orange	INTG156403	<b>\$819.00</b>	<b>\$475.00</b>
PIPETBOY Pro Blue	INTG156404	<b>\$819.00</b>	<b>\$475.00</b>

*Pipettor complete with spare sterile filter, two batteries, mains adaptor and wall mount.*

**Only  
\$475  
Each!**



Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

**VALID UNTIL:  
30 September 2019**

**Ref. Code:  
INTG-PRO**

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.