

PHCbi CO2 Incubator Package

Optimising cell culture outcomes and reproductivity.

PHCbi Incusafe CO2 Incubators are the result of more than 50 years of technological innovation.

Features:

- UV sterilisation of circulating air and water pan
- H2O2 decontamination (maximum 3-hour process)
- Dual IR CO2 sensor – no need for removal during decontamination
- Copper embedded stainless-steel interior
- Integrated shelf supports for easy clean up

Scope of Package: 2x PHCbi MCO-170AICUVHL-PE CO2 incubators, 1x MCO-170PS stacking kit, 1x roller base and also includes 1x H2O2 generator + 6 pack of H2O2 reagent for H2O2 decontamination.

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
2x SANYMCO-170AICUVHL-PE 1x SANYMCO-170PS 1x SANYMCO-170AICTROL 1x SANYMCO-HP 1x SANYMCO-H2O2	SANYCO2PACKAGE1	\$18,644.00	\$12,990.00

Items are available ex stock at time of print.

**\$12,990
for Stacked
CO2 Incubator
Package**



Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

VALID UNTIL:
30 September 2019

Ref. Code:
SANY-CO2

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.