

New PHCbi Pharmaceutical Refrigerators from Bio-Strategy

Delivering the same stable and reliable temperature control, in a new look.

Top mounted air circulation blowers offer fast recovery times following door openings, even with frequent access and high ambient temperatures. The new temperature display will display the interior temperature in 0.1°C increments.

Two size options available with either standardised wired shelving or ergonomic draws. There's a unit available to meet any specific need.

The PHCbi pharmaceutical refrigerator range features:

- Forced cool air circulation assures quick recovery and temperature uniformity
- Precise temperature control
- Cycle defrost function with no irregular temperature spikes
- 1 or 2, double pane glass doors (self closing)
- New optional black out panels

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

Description	Code
Refrigerator Pharmaceutical 684L	SANYMDF-MPR-722
Refrigerator Pharmaceutical 671L Full Roller Drawers	SANYMDF-MPR-722R
Refrigerator Pharmaceutical 1,364L	SANYMDF-MPR-1412
Refrigerator Pharmaceutical 1,359L Full Roller Drawers	SANYMDF-MPR-1412R

**Available
Q4 2019!**



Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.