

# Scilabware MBL® Volumetric Flasks

## Volumetric Class A flasks – buy now and save!

**Class A volumetric flasks printed with blue enamel high contrast markings.**

- Complies with ISO 1042 and DIN 12664
- Finely calibrated to Class A tolerances
- Each flask is laser etched with an individual serial number, batch code and date of manufacture
- Batch certificate available to download online
- Manufactured from chemically resistant borosilicate glass
- Conformity 'H' certified
- Supplied with polyethylene stopper

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)**

Capacity (ml)	Tol. (+/- ml)	Approx Height (mm)	Stopper Size	Pack Size	Availability	Code	Sale Price (excl. GST)
<b>Flasks Volumetric Class A with PE Stopper</b>							
5	0.025	70	10/19	5/Pkt	In Stock	BIBBFRD164	<b>\$29.34</b>
10	0.025	90	10/19	5/Pkt	In Stock	BIBBFRD166	<b>\$24.26</b>
20	0.04	110	10/19	5/Pkt	In Stock	BIBBFRD167	<b>\$25.06</b>
25	0.04	110	10/19	5/Pkt	In Stock	BIBBFRD168	<b>\$25.09</b>
50	0.06	140	12/21	5/Pkt	In Stock	BIBBFRD170	<b>\$30.19</b>
100	0.10	170	14/23	5/Pkt	In Stock	BIBBFRD172	<b>\$34.93</b>
200	0.15	210	14/23	5/Pkt	Enquire	BIBBFRD174	<b>\$53.36</b>
250	0.15	220	14/23	2/Pkt	Enquire	BIBBFRD176	<b>\$24.45</b>
500	0.25	260	19/26	2/Pkt	Enquire	BIBBFRD178	<b>\$28.32</b>
1000	0.40	300	24/29	2/Pkt	Enquire	BIBBFRD180	<b>\$40.53</b>
2000	0.60	370	29/32	1/Pkt	Enquire	BIBBFRD182	<b>\$40.33</b>
5000	1.20	475	34/35	1/Pkt	Enquire	BIBBFRD184	<b>\$204.42</b>

*If not In Stock delivery approximately 6-8 weeks, unless already on order.*

**VALID UNTIL:  
30 September 2019**

**Ref. Code:  
BIBB-FLASK**

**OVER  
25%  
OFF!**



Bio-Strategy Pty Ltd.

T: 1800 008 453  
E: [sales.au@bio-strategy.com](mailto:sales.au@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
<https://shop.bio-strategy.com>

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.