

Tomy SX-700E Top Loading Electrically Heated Autoclave

Includes one free preventative maintenance inspection!

Compact Mobile Floor Standing Autoclave that does not require plumbing!

The Tomy SX-700E top loading autoclave has a 79-litre capacity and is supplied with two baskets. Three operating temperature ranges: warming 45 to 95°C, heating 45 to 104°C and sterilising 105 to 135°C. Safety devices: water level sensor, current leakage breaker, lid interlock, overheat prevention, over pressure prevention, open temperature sensor detection and safety valve. With integrated (but removable) condensate bottle to trap condensate from chamber during venting.

- Simple lid opening using either hand or foot
- Cooling fan to reduce run time
- One Free Preventative Maintenance visit (metro area's only)
- Supplied with two baskets
- Digital temperature and time display with analogue gauge for pressure
- Five sterilisation profiles to select from (in addition to normal sterilisation)

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 1800 008 453 or email sales.au@bio-strategy.com

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Autoclave 79L*	TOMYTMYSX-700E	\$13,725.00	\$9,500.00

*Requires 240V 50Hz 15-amp supply.

On Sale for only \$9,500!



Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
<https://shop.bio-strategy.com>

**VALID UNTIL:
30 September 2019**

**Ref. Code:
TOMYSX-700E**

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.