

## New Range of UVP Penray Lamps!

VALID UNTIL:  
30 September 2020


Ref. Code:  
UPI-PENRAY

Invented by UVP, now Analytik Jena, the Pen-Ray Lamps are manufactured from the highest grade of quartz material. The design places the illumination on one end of the lamp, making it ideal to use as a reference source for small apertures or in calibration applications for instruments.

- Available in standard sizes and wavelengths
- Customisation options are available (quartz type and size)
- Power supply not included, must be ordered separately

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
Model 11SC-1. 254nm, light length = 53.8mm, Ø 6.5mm. One of the most popular Pen-Ray Lamps	UVPI90-0012-01	\$595	\$500
Model 11SC-1L 365nm, light length = 57.2mm, Ø 6.5mm	UVPI90-0019-01	\$640	\$535
Model 11SC-2 365nm, light length = 19.1 mm, Ø 6.5mm	UVPI90-0020-01	\$570	\$480
Model 3SC-9 254nm, light length = 228.6 mm, Ø 9.5mm	UVPI90-0004-01	\$1,155	\$970
Recommended Power Supply (11SC-1/11SC-1L/11SC-2)	UVPI99-0055-02	\$685	\$575
Recommended Power Supply (3SC-9)	UVPI99-0004-02	\$675	\$565

**Perfect for  
calibration  
applications!**



Bio-Strategy Pty Ltd.

T: 1800 008 453  
sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com