

# New Range of UVP UV Bench Lamps!

VALID UNTIL:  
30 September 2020

Ref. Code:  
UVP-BENCH

UVP XX-Series UV Bench Lamps - various fixture sizes, watts and wavelengths are available. Lamps include holes on each end of the fixture for hanging or mounting the lamp. Lamps are manufactured with a durable aluminum housing and powder coated scratch resistant finish.

Shortwave lamps can be used in an upward facing position for purifying air or in a downward position for sanitizing work surfaces 254 nm UV can kill up to 99% of organisms from bacteria, yeast, mold spores, germs and viruses.

- Select from 15 watt, 20 watt or 40 watt options
- Select from 254 nm, 365 nm, 365 nm BLB and 302 nm.
- XX-15 lamp dimensions LxWxH 502 x 152 x 108 mm, 2.3 kg
- XX-20 lamp dimensions LxWxH 610 x 152 x 108 mm, 4.5 kg
- XX-40 lamp dimensions LxWxH 1,260 x 152 x 108 mm, 9.0 kg



**Protect against  
contamination!**

Description	Code	Price (excl. GST)	Sale Price (excl. GST)
XX-15S 15W Bench Lamp 254nm	UVPI95-0042-09	\$800	\$670
XX-15M 15W Bench Lamp 302nm	UVPI95-0042-15	\$840	\$705
XX-15L 15W Bench Lamp 365nm	UVPI95-0042-13	\$775	\$650
XX-15BLB 15W Bench Lamp 365nm BLB	UVPI95-0042-11	\$800	\$670
XX-20S 20W Bench Lamp 254nm	UVPI95-0045-08	\$965	\$805
XX-20BLB 20W Bench Lamp 365nm BLB	UVPI95-0045-05	\$965	\$805
XX-40S 20W Bench Lamp 254nm	UVPI95-0043-08	\$1,095	\$915
XX-40M 20W Bench Lamp 302nm	UVPI95-0043-14	\$1,065	\$890
XX-40BLB 20W Bench Lamp 365nm BLB	UVPI95-0043-05	\$1,065	\$890
XX-15 Exposure Stand	UVPI18-0062-01	\$915	\$770
Filter for XX-15S	UVPI38-0172-02	\$1,535	\$1,290

Bio-Strategy Pty Ltd.

T: 1800 008 453  
sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com