

Duran Laboratory Bottles on Special!

VALID UNTIL:
30 September 2020

Ref. Code:
SCOT-LAB

Duran Laboratory Bottles (aka Schott) are manufactured in Duran 3.3 borosilicate glass. Each bottle has its own 8 character retrace number with certificate downloadable online.

- Easy to read, white ceramic fired on durable markings
- Glass Type I conforms to USP/EP standards
- High temperature, thermal shock resistant and uniform wall thickness
- Blue PP autoclavable screw cap and pouring ring, for drip free, safe work +140°C

**Great
Value!**



Description	Pack size	Code	Price	Sale Price
BOTTLE LAB DURAN 50ML GL32 BLUE CAP AND RING	1 item	SCOT218011753	\$9.75	\$8.00
BOTTLE LAB DURAN 100ML GL45 BLUE CAP AND RING	1 item	SCOT218012458	\$8.80	\$6.90
BOTTLE LAB DURAN 250ML GL45 BLUE CAP AND RING	1 item	SCOT218013651	\$9.98	\$8.70
BOTTLE LAB DURAN 500ML GL45 NO CAP AND RING	1 item	SCOT218014401	\$9.78	\$7.50
BOTTLE LAB DURAN 500ML GL45 BLUE CAP AND RING	1 item	SCOT218014459	\$11.17	\$8.70
BOTTLE LAB DURAN 1L GL45 BLUE CAP AND RING	1 item	SCOT218015455	\$15.70	\$13.70
BOTTLE LAB DURAN 2L GL45 BLUE CAP AND RING	1 item	SCOT218016357	\$32.59	\$28.50
BOTTLE LAB DURAN 5L GL45 BLUE CAP AND RING	1 item	SCOT218017353	\$203.62	\$178.00
BOTTLE LAB DURAN 20L GL45 BLUE CAP AND RING	1 item	SCOT218019157	\$720.22	\$585.00

Bio-Strategy Pty Ltd.

T: 1800 008 453
sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com