

Kimble Tubes on Special!

VALID UNTIL:
30 September 2020

Ref. Code:
KIML-TUBES

These premium quality tubes from Kimble all have sturdy uniform bottoms and consistent lengths.

- Borosilicate glass tubes are manufactured using 51 expansion borosilicate glass and conform to USP Type I and ASTM E438 Type I, Class B requirements. They offer excellent chemical resistance against sodium leaching.
- Soda-Lime tubes are manufactured using 90 expansion soda-lime glass and conform to USP Type III requirements.

**From \$9.50
per pack!**



Description	Pack size	Code	Price	Sale Price
TUBE PLAIN DISPOSABLE BOROSILICATE 1ML 6X50MM (DURHAM)	100 items	KIML73500-650	\$16.15	\$9.50
TUBE PLAIN DISPOSABLE BOROSILICATE 10ML 13X100MM	1000 items	KIML73500-13100	\$115.00	\$90.00
TUBE PLAIN DISPOSABLE BOROSILICATE 15ML 16X100MM	1000 items	KIML73500-16100	\$170.00	\$134.00
TUBE PLAIN DISPOSABLE BOROSILICATE 19ML 16X125MM	1000 items	KIML73500-16125	\$150.00	\$129.00
TUBE PLAIN DISPOSABLE BOROSILICATE 23ML 16X150MM	1000 items	KIML73500-16150	\$129.00	\$85.00
TUBE PLAIN DISPOSABLE BOROSILICATE 28.5ML 18X150MM	500 items	KIML73500-18150	\$125.00	\$85.00
TUBE PLAIN DISPOSABLE SODA LIME 23ML 16X150MM	1000 items	KIML60GM166	\$139.00	\$94.00

Bio-Strategy Pty Ltd.

T: 1800 008 453
sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com