

Cytiva kits for use in COVID-19 testing!

**Contact us
for a trial!**

VALID UNTIL:
30 September 2020

Ref. Code:
CYTIVA-KITS

RNAspin Mini Kit:

Designed for rapid extraction of high-quality RNA from a wide range of sample types. Includes prefilters and DNase I.

RNAspin Mini RNA Isolation Kit is a complete RNA purification kit.

- High-quality output RNA from diverse sample types
- Recovers high-quality total RNA since genomic DNA is removed via on-column DNase I treatment
- For maximum yield and purity, prefilters are included to reduce lysate viscosity
- Results can be obtained with even small amounts of sample
- Lysis buffer is less susceptible to foaming to ensure valuable RNA sample is not wasted
- Simple and convenient format suitable for all levels of expertise



Illustra Ready-To-Go RT-PCR Beads:

Pre-dispensed, freeze-dried beads that include the reagents necessary for one-step reverse transcription-PCR with high sensitivity and reproducibility.

- Preformulated, predispensed, single-dose, ambient-temperature-stable beads ensure greater reproducibility between reactions, minimize pipetting steps, and reduce the potential for pipetting errors and contamination.
- Optimized as a one-tube, one-step RT-PCR reaction for both cDNA synthesis and PCR.
- Each lot of RT-PCR Beads is function tested for its ability to generate highly specific PCR products to ensure lot-to-lot reproducibility.
- For highly reproducible one-step reverse transcription-PCR.
- Ready-To-Go RT-PCR Beads are stable at room temperature and designed for performing single-tube one-step reverse transcription-PCR. Each room-temperature-stable bead contains M-MuLV Reverse Transcriptase, RNase Inhibitor, buffer, nucleotides, and Taq DNA Polymerase. The only additional reagents required are water, template RNA, and primers. The reagents are optimized for full-length cDNA synthesis to >7.5 kb and optimal sensitivity from PCR.



cytiva
Distributor



Description	Code	Price	Sale Price
illustra RNAspin Mini Kit	25-0500-87	\$91.85	Free for Trial
illustra RNAspin Mini Kit,	25-0500-70	\$140.85	\$130.95
illustra RNAspin Mini Kit,	25-0500-71	\$293.00	\$272.50
illustra RNAspin Mini Kit	25-0500-72	\$1,261.00	\$1,172.70
illustra RNAspin 96 Kit	25-0500-75	\$1,415.00	\$1,315.95
illustra Ready-To-Go RT-PCR Beads (0.2 ml Hinged Tube with Cap)	27-9259-01	\$516.00	\$479.90
illustra Ready-To-Go RT-PCR Beads (0.5 ml Tubes)	27-9266-01	\$541.00	\$503.10

Bio-Strategy Pty Ltd.

T: 1800 008 453
sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com