

Great Value Packages from Labconco!

VALID UNTIL:
30 September 2020

Ref. Code:
LBCP-PACK

Choose Labconco for preservation and concentration with the FreeZone Freeze Dryers and Centrivap Centrifugal Vacuum Concentrators.

FreeZone Freeze Dryers (lyophilizers) are designed to handle the lyophilization needs of research and pilot plant laboratories. Select from benchtop or mobile console options. A 4.5 litre lyophilizer that reaches -105°C is also available, suitable for samples containing dilute methanol and ethanol. Collectors that reach -84°C are ideal for freeze drying samples with low eutectic temperatures including ones containing acetonitrile.

- Collector options: -50°C, -84°C or -105°C
- Capacities ranging from 2.5 L to 18 L
- Stainless steel or PTFE-coated chambers and accessories
- Lyo-Works™ OS provides a real time display of collector temperature, vacuum level and optional End-Zone™ end point detection system.

Centrivap Centrifugal Vacuum Concentrators use a combination of centrifugal force, vacuum and heat to speed laboratory evaporation of multiple small samples. CentriVap Cold Traps protect your vacuum pump from corrosive chemicals.

- Optional 300 W heat boost which surrounds the chamber wall, to further speed evaporation of high boiling point solvents
- Select cold traps at -50°C, -84°C or -105°C, with optional traps downstream to help protect your pump for corrosive chemicals
- Stainless steel or PTFE-coated chambers and accessories are available

Let Bio-Strategy
configure you a
system for your
application!



Description	Code	Price	Sale Price
1x LBCP700201040 FreeZone 2.5L -50° C Benchtop Freeze Dryer with stainless steel coil and collector 1x LBCP7201700 4-Port Manifold, 8.6" w x 8.6" d x 11.2" h 1x LBCP7739402 Vacuum Pump with Filter Exhaust, 117 liters/min., with IEC plug 2x LBCP7540800 600 ml Complete Fast-Freeze Flask 2x LBCP7541100 750 ml Complete Fast-Freeze Flask 4x LBCP7547600 Fast-Freeze 45° Adapter (Stainless Steel)	PACKAGE1	\$26,513	\$20,995
LBCP7810040 CentriVap Benchtop Concentrator with Glass Lid LBCP7811038 CentriVap -50 Cold Trap with stainless steel lid LBCP7450700 CentriVap Hexagonal Rotor, 1.5 ml tubes LBCP7739402 Vacuum Pump with Filter Exhaust, 117 liters/min.,with IEC plug	PACKAGE2	\$28,875	\$22,390

Bio-Strategy Pty Ltd.

T: 1800 008 453
sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com