

Save on Miele High Quality Washer/Disinfector Package

VALID UNTIL:
30 September 2020

Ref. Code:
MIEL-PG8504

The PG8504 glassware washer incorporates numerous features to guarantee optimum results for every wash. The washer boasts an improved capacity, allowing you to clean large number of utensils and glassware quickly. Freestanding or under counter, the stainless steel PG8504 has been specially constructed to reduce heat and noise emissions and the interior has been welded together using lasers. The PG8504 is designed with convenience in mind, utilising simple to load door dispensers. The door also integrates the TouchControl display, allowing for a quick and simple selection of 6 different programs, and is built in to the handle for ease.

Please note:

Customer to organise their local plumber and electrician to connect unit to services

- Door dispensers for powder detergent and neutraliser
- Powerful internal wash pump (500L/min) with integrated heating elements ensures a thorough clean while reducing the chance of contamination
- Free-standing or under-bench design: fits in space 600mmW x 600mmD x 850mmH (820mm with lid removed)
- The 3-line display keeps the operator fully informed of the washer status
- 6 programs and temperature range of up to 70°C
- Suitable for operation on 240V 50Hz 13 amp supply



Description	Code	Price	Sale Price
PG8504 Package 1 Washer with upper + lower baskets and accessories	PG8504-Package 1	\$10,893	\$9,995

**Save
over
\$800!**

Bio-Strategy Pty Ltd.

T: 1800 008 453
sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com