

# Configurable Autoclaves by Priorclave

Priorclave have manufactured autoclaves in London for over 30 years. Built to rigid British standards – the units comply to AS/NZS standards and the RSC (rectangular section chamber) models feature capacities from 230 to 700 litres.

The units can be configured to suit a multitude of applications with options including manual or power doors, water/steam jackets, load sensed process timing (wandering probes), pre-cycle vacuum and vacuum cooling, vacuum drying, trolleys, loading cassettes, effluent and vent filters, condensers, printers, containers and trivets, air ballasting etc .

- Capacities: 230, 285,350, 400, 450 and 700 litres.
- Front Door and Pass Through models to suit your requirements.
- 10 year chamber and door(s) warranty (manufactured from 316SS).
- One flexible program with 5 and 10 program options.
- USB port, media warming facility, timed free-steaming, external fan cooling and two shelves supplied as standard.
- Elements in chamber/internal steam generator or external steam generator options.



Contact us to  
discuss your  
Autoclave  
requirements!

Bio-Strategy Pty Ltd.

T: 1800 008 453  
sales.au@bio-strategy.com  
www.bio-strategy.com  
shop.bio-strategy.com