

TOMY Autoclave on Special

VALID UNTIL:
30 September 2020

Ref. Code:
TOMY-SX700E

The Toploading SX-700E Autoclave from Tomy (Japan) requires no costly plumbing or electrical services as the unit has a self contained reservoir for exhaust steam and can be run off a 240V 50Hz 15 amp supply.

The microprocessor controller with digital display keeps the end users informed of the current instrument status. The foot pedal mechanism assists with the lid opening and the unit has safety mechanisms including lid interlock and water level detection. Supplied complete with 3 x ACA-700B mesh baskets.

- 2-year parts and labour warranty
- Free preventative maintenance visit at the end of the first year to ensure unit is meeting manufacturer's specifications.



**Includes
preventative
maintenance!**

Description	Code	Price	Sale Price
Tomy Model SX-700E Toploading 79 litre Autoclave with 2 baskets.	TMYSX-700E	\$17,086	\$12,995
Additional Stainless Steel Basket for SX-700E Autoclave	132529027	\$688	included

Bio-Strategy Pty Ltd.

T: 1800 008 453
sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com