

Heidolph Vortex Mixers

Reax Top / Reax Control Vortex Mixer

- The standard or precision model
- The shaking orbit of 5mm reliably and quickly achieves an even distribution
- Reax control with electronic speed control – the speed remains constant even in the low range and at load changes

Multi Reax

- Shaking up to 26 samples simultaneously and achieving excellent mixing results
- Specific adaptors are not required
- Excellent mixing results are achieved with the 3mm shaking orbit, even with large samples with solid content
- Speed adjustment from 150 to 2,000rpm
- Timer function up to 999 minutes

Incredibly strong and robust, the strong shaking motion guarantees excellent mixing results without exception.

- Designed and made in Germany: 3-year warranty on all devices and an average operational lifespan of 10-years
- A worthwhile investment: maintenance and spark-free motors and are excellently suited for years of continuous use
- The sealed housing reliably protects against corrosion
- Fast, even distribution, even with solids and highly viscous media: ideally suited for short- or long-term operation



VALID UNTIL:
31 December 2020

Ref. Code:
HEID-MIX



In Stock!
3-Year Warranty!
Lifespan of 10-Years!



**MADE IN
GERMANY**

3-year warranty on all devices
and an average operational
lifespan of 10 years

Multi-stage quality checks in
development and production

Premium service according to the
"research made easy" principle

Description	Code	Price (ex GST)	Sale Price (ex GST)
Reax Top with Small Pick-Up Attachment	HEID541-10000-00	\$563.00	\$450.00
Reax Control with Small Pick-Up Attachment	HEID541-11000-00	\$769.00	\$615.00
Multi Reax	HEID545-10000-00	\$2,612.00	\$2,089.00

Bio-Strategy Pty Ltd.

T: 1800 008 453
sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com