

Kimble Clearance Sale

Super Special Prices!

While stocks last.

KIMBLE® Plain Disposable Borosilicate Glass Tubes:

- 51 expansion borosilicate glass offers excellent chemical resistance against sodium leaching, a factor common in soda lime glass products
- Tubes are packed in convenient trays for ease of use on bench tops and in drawers
- The 6 x 50 mm size is often referred to as a “Durham tube”
- Trays are shrink-wrapped in plastic film, with tubes remaining in constant alignment for protection against in-transit breakage and contamination
- No marking spot
- Designed from ASTM Specification E890
- Manufactured from 51 expansion borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class B requirements

ON SALE!



Description	Pack Size	Code	Sale Price (excl. GST)
Tube Plain Durham 1ml Disposable Borosilicate 6x50mm	100 Items	KIML73500-650	\$5
Tube Plain Disposable Borosilicate 23ml 16x150mm	1,000 Items	KIML73500-16150	\$50
Tube Plain Disposable Borosilicate 28.5ml 18x150mm	500 Items	KIML73500-18150	\$50

Contact us for more information
and to place your orders

Bio-Strategy Ltd.

T: 0800 34 24 66
E: sales.nz@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2021

Ref. Code:
KIML-SALE

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.