

Heidolph Hei-Vap Precision Package

Solution for the most demanding evaporation tasks!

A motor-lift model with built-in vacuum control and AutoAccurate sensor automatically detects the boiling point temperature of the solvent and adjust the vacuum.

The Hei-VAP Precision provides you with the maximum range of parameters, including but not limited to, reproducible applications and fully automated processes. Even solutions with different boiling points can be identified and distilled.



Price valid until
stock lasts!

| Description | Code | Price (excl. GST) | Sale Price (excl. GST) |
|--|--------------------|----------------------|----------------------------|
| Hei-Vap Precision ML G3B With Built-In Vacuum Control, Motor Lift, Coated Vertical Glassware. 1L Evaporation Flask and 1L Receiving Flask are Included | HEID564-01310-00 | \$10,196 | included in package |
| Rotovac Vario Control Vacuum Pump | HEID591-00141-00 | \$8,789 | included in package |
| Condenser for Rotovac Vario Control Vacuum Pump | HEID591-00084-00 | \$1,094 | included in package |
| Temperature AutoAccurate Sensor | HEID569-00040-00 | \$449 | included in package |
| Tube Set | HEID591-35000-00 | \$406 | included in package |
| Hei-Vap Precision ML G3B Package (Package Includes: Hei-Vap Precision ML G3B with Built-In Vacuum Control, Motor Lift, Coated Vertical Glassware; Ex-Demo Unit, Warranty: 2-Years) | HEIDHEI-VAPPACKAGE | \$20,934 | \$13,693 |

Heidolph Evaporators:

- Highest safety standards
- 10-year lifespan
- Vapour tube with patented clamping sleeve: Heidolph patented clamping sleeve system allows easy removal of the vapour tube from the drive and it will reduce your spare parts budget up to 75%
- Patented flask clamp Easy-Clip: remove evaporation flask using one hand, no glass-to-glass contact
- Additional unique features designed by experts in evaporation
- Ideal for foaming samples
- Multiple solvents samples
- Solvent recycling
- Toxic samples
- Large 4.3" colour graphic display which features all parameters, integrated vacuum controller and includes programs for automated distillations
- Save time with your individualised pre-programmed parameter settings for your 30 most common applications used on a daily basis by saving them as "Favourites"
- Each ramp offers 20 separate steps for demanding applications
- Apply an individual name to those ramps to avoid a mix-up
- Experience the advantages of the automatic process timer which turns off your evaporator at a pre-programmed time by raising the flask on motor lift models and releasing vacuum
- Hei-VAP Precision models feature a built-in vacuum controller as a standard
- For increased safety all bath temperature adjustments can be limited
- Even from a distance you will be able to monitor the numbers on the display easily

Intelligent evaporation with Heidolph rotary evaporators: overheat protection design that prevents water baths running dry, provides a detachable panel for superior ease of use from outside closed fume hoods and offers a proven vacuum seal and vapour tube system guaranteed for years of operation!

Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453

E: sales.au@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com



[linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2021

Ref. Code:
HEID-VAP

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.