

AcroPrep™ 24-Well Filter Plates

For protein purification and 0.2µm sterile filtration workflows.

New Innovation in a 24-well Format for Workflow Improvements

Pall Laboratory has two 24-well, 7ml, filter plates for use with protein purification and general sterile filtration workflows. Utilising Pall's proprietary high-performance multi-layer filter media and membranes, these plates offer time savings, strong performance claims and streamlined workflow improvements in a 24-well plate format. Each filter plate comes with a V-shaped bottom collection plate and lid.

Cell Clarification and Sterile Filtration Plate for Protein Purification Workflows

Pall's dual media clarification and sterile filtration AcroPrep 24-well filter plate can both clarify and 0.2µm sterile filter in a single device and workflow step. The top layer features Pall's Seitz® depth media and efficiently captures whole cells and cell debris. The lower Supor® EKV layer provides an efficient sterile filtration layer. The combination effortlessly recovers proteins from whole cell cultures with varying viabilities of up to 25M+ cells/ml in minutes using either a vacuum manifold or plate centrifuge.

0.2µm Filtration Plate for Sterile Filtration Workflows

Pall's sterile filtration AcroPrep 24-well filter plate is well suited for plate-based sterile filtration needs such as for media, reagent, serum or proteins. Compatible with either a vacuum manifold or plate centrifuge, the filter plate utilises Pall's Supor® EKV membrane which is a dual-layer PES 0.65µm membrane integrated with a highly asymmetric 0.2µm membrane for fast, efficient, and sterilising grade filtration. The plate is tested to meet sterile *Brevundimonas diminuta* (Bdim) retention claims.



Description	Pack Size	Code
AcroPrep 24W Clarification Filtration Plate Sterile	2 Items	PALL97016
AcroPrep 24W Sterile Filtration Plate	2 Items	PALL97017
AcroPrep 24W Clarification Filtration Plate Sterile	8 Items	PALL97026
AcroPrep 24W Sterile Filtration Plate	8 Items	PALL97027

NEW!



Contact us for more information
and to place your orders

Bio-Strategy Pty Ltd.

T: 1800 008 453
E: sales.au@bio-strategy.com
www.bio-strategy.com
shop.bio-strategy.com

 [linkedin.com/company/bio-strategy](https://www.linkedin.com/company/bio-strategy)

VALID UNTIL:
31 March 2021

Ref. Code:
PALL-ACRO

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.